

Cristobal⁺
AND
Calibra[®]
esthetic resin cement

Chairside Preparation and Bonding Guide

The following preparation and bonding guidelines for restorative indications of this system will ensure the clinical success of your restorations.

Indications

- Metal-free/fiber-free
 - Inlays/Onlays:
All teeth
 - Veneers:
All teeth
 - Crowns:
All teeth through 1st molar
 - Bridges:
All teeth through 2nd bicuspid,
one pontic bridge only
- With metal substructure
 - Crowns & Bridges
(w/metal substructure):
All teeth
- Implants
 - All teeth/bonded-to-metal

Contraindications

- Do not use with fiber substructures


Preparation Guidelines

General specifications:

1. Use a chamfer margin around the gingival margins. Margins that extend to the occlusal or incisal must be a 90° or butt joint.
2. Avoid sharp angles.
3. Maintain an even reduction of anatomical form.
4. These are minimum reductions. Strength and opacity will be increased with more tooth reduction.

Note: Modifications to the above guidelines must be made with caution, using clinical judgment.

Inlay/Onlay Preparation Protocol




Veneer Preparation Protocol



* 0.5 is the minimum tooth reduction. Added strength and opacity may be increased by greater tooth reduction.

Metal-Free Preparation Protocol



NOTE: Metal-free bridge proximal connector area 1.5 to 2.0-mm reduction

Composite-Bonded-To-Metal Preparation Protocol



See complete directions for use for indications and contraindications.

Bonding

Cristobal+ restorations require a simple clinical protocol. Providing the technician with the shade of the prepared tooth as well as the final shade will ensure the desired result and simplify cementation.

Selecting Calibra® Esthetic Resin Cement and Prime & Bond® NT™ Dual Cure Dental Adhesive will help ensure your success with Cristobal+ restorations.

Always follow the manufacturers' directions.

An abbreviated technique follows.

Step 1: Remove the provisional restoration. Clean and dry the site.

Step 2: Verify the fit of the restoration utilizing a try-in paste of glycerin.

Step 3: Isolate the prepared tooth with a rubber dam and retainer.

Step 4: Apply and light-cure Prime & Bond® NT™.

Step 5: Treat the internal aspects of the restoration with a silane coupling agent.

Step 6: Apply Calibra® Esthetic Resin Cement to the surface of the restoration.



Step 7: Seat the restoration and remove gross excess cement.

Step 8: Light-cure the entire restoration from all surfaces.

Step 9: Adjust and finish with fine diamonds and polishers.

For more information on the Calibra® Esthetic Resin Cement, contact your local dental distributor or Dentsply Caulk at 800-532-2855.

Calibra Complete Kit Item # 607060
Calibra Operatory Kit Item # 607059



Cristobal+ AND **Calibra**
esthetic resin cement

DENTSPLY
CERAMCO

ISO 9001 CERTIFIED

Six Terri Lane • Burlington, NJ 08016 • 800-487-0100 • In Canada: 1-888-881-0886 • www.ceramco.com • ©2000 Ceramco Inc. All Rights Reserved.